

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2-(Dimethylamino)ethyl methacrylate - contains 700-1000 ppm monomethyl ether hydroquinone as inhibitor, 98%

Product Number:

234907

CAS Number:

2867-47-2

MDL:

MFCD00008589

Formula:


C₈H₁₅NO₂

Formula Weight:

157.21 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Proton NMR Spectrum

Conforms to Structure

Purity (GC)

≥ 97.5 %

MEHQ as Inhibitor

700 - 1000 ppm

Specification: PRD.1.ZQ5.10000018208

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.