

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Monoclonal Anti-Cytokeratin, pan antibody produced in mouse - clone PCK-26, ascites fluid

Product Number:

C1801

**TEST****Specification**

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Very Faint Yellow to Very Dark Yellow to Very Dark Brown
Protein by Biuret	10 - 55 mg/ml
Mouse IgG1 by RID	7.0 - 14.0 mg/ml
IHC - Tissue	Human placenta or skin
IHC - Fixation	FFPE-Formalin Fixed Paraffin Embedded
IHC - Titer	1:300-1:600
Recommended Retest Period	-----
5 Years	

Specification: PRE.0.ZQ5.10000015219

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.