

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Sodium carbonate – anhydrous, powder, 99.999% trace metals basis

Product Number:

451614

CAS Number:

497-19-8

MDL:


MFCD00003494

Formula:

CNa<sub>2</sub>O<sub>3</sub>

Formula Weight:

105.99 g/mol



| TEST                                   | Specification      |
|----------------------------------------|--------------------|
| Appearance (Color)                     | White              |
| Appearance (Form)                      | Powder             |
| ICP Major Analysis                     | Confirmed          |
| Confirms Sodium Component              |                    |
| Purity                                 | Meets Requirements |
| 99.999% Based On Trace Metals Analysis |                    |
| Trace Metal Analysis                   | ≤ 15.0 ppm         |

Specification: PRD.0.ZQ5.10000068134

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.