


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
18 α -Glycyrrhetic acid - $\geq 95\%$

Product Number: G8503
CAS Number: 1449-05-4
MDL: MFCD00064897
Formula: C₃₀H₄₆O₄
Formula Weight: 470.68 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
EmM in EtOH	10.9 - 11.6
Wavelength	243 - 245 nm
Purity (TLC)	$\geq 95\%$

Specification: PRD.1.ZQ5.10000026239

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.