

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification


Product Name:
Glyoxal solution - 40 wt. % in H₂O

Product Number: 128465

MDL: MFCD00006957

Formula: C₂H₂O₂

Formula Weight: 58.04 g/mol



TEST

Specification

Appearance (Color)

Material May Darken In Storage With No Loss Of Purity

Colorless to Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Titration by Hydroxylamine

38.0 - 42.0 %

Specification Date : 01/31/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.