

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Fusel oil - natural, FGProduct Number: **W249715**

CAS Number: 8013-75-0

Formula: C₅H₁₂O

TEST	Specification
------	---------------

Appearance (Color)	Colorless to Light Brown and Colorless to Light Orange and
Colorless to Light Red	
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.385 - 1.410
Infrared spectrum	Complete
Specific Gravity at 20°C	0.785 - 0.860
Arsenic (As)	≤ 3 ppm
Mercury (Hg)	≤ 1 ppm
Cadmium (Cd)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period 60 Months	-----

Specification: PRD.1.ZQ5.10000027753

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.