

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

8-Chloroadenosine 3',5'-cyclic-monophosphate - ≥85%

Product Number:

C1081

CAS Number:

41941-56-4

MDL:

MFCD00153938

Formula:


C₁₀H₁₁ClN₅O₆P

Formula Weight:

363.65 g/mol

Storage Temperature:

-20 °C



TEST

Specification

Appearance (Color)

White to Off-White

Appearance (Form)

Powder

% Purity (HPLC)

≥ 85.0

% Water

≤ 13.5

Solubility (Color)

Colorless to Yellow

Solubility (Turbidity)

Clear

1 mg + 0.1 mL of 1M NH₄OH

Specification: PRD.2.ZQ5.10000054761

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.