


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
DCA Deblock (0.36M dichloroacetic acid in toluene)

Product Number: 707066



Formula: C<sub>2</sub>H<sub>2</sub>Cl<sub>2</sub>O<sub>2</sub>  
Formula Weight: 128.94 g/mol

TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum ATR	Conforms to Structure
Color Test	≤ 20 APHA
Purity Titration by NaOH (M)	0.34 - 0.38 M
Water (by Karl Fischer)	≤ 100 ppm
Turbidity	≤ 3.5 NTU

Specification: PRD.0.ZQ5.10000027126

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.