SIGMA-ALDRICH®


sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name: 1,4-Dibenzoylbenzene - 99%

Product Number: CAS Number: MDL: Formula: Formula Weight: 445932 3016-97-5 MFCD00014087 C20H14O2 286.32 g/mol



TEST

Appearance (Color) Appearance (Form) Infrared spectrum Purity (HPLC)

Specification: PRD.1.ZQ5.1000005161

Specification

White to Light Yellow to Light Brown Powder Conforms to Structure > 98.5 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.