


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
H-Gly-His-Lys-OMe trifluoroacetic acid salt - 95%

Product Number: 776629



Formula: C₂₁H₂₉F₉N₆O₁₀
Formula Weight: 696.47 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	White to Faint Yellow and Faint Beige and Faint Brown
Appearance (Form)	Powder or Crystals
Proton NMR spectra	Conforms to Structure
Optical Rotation c = 0.5%, MeOH	-11.0 - -5.0 deg
Carbon	34.2 - 38.2 %
Nitrogen	11.4 - 12.7 %

Specification: PRD.1.ZQ5.10000048326

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.