

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Silver nanoplates – ~650 nm (Resonant), ~1 mg/mL (in water with 5 mM sodium borate buffer), PVP functionalized
Product Number: 807524

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	Conforms to Requirements
Appearance (Form)	Liquid
Wavelength	600 - 699 nm
ICP Major Analysis	Confirmed
Confirms Silver Component	
Concentration	
Mass(Ag)	
Specification: Report Only	
Diameter	
(TEM)	
Specification: Report Only	
Standard Deviation	
Specification: Report Only	
pH	
Specification: Report Only	
Particle Surface	Confirmed
PVP	
Recommended Retest Period	-----
1 Year	
Specification: PRD.0.ZQ5.10000073577	