


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
trans-Cinnamic acid - ≥99%

Product Number: C80857
CAS Number: 140-10-3
MDL: MFCD00004369
Formula: C₉H₈O₂
Formula Weight: 148.16 g/mol



TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Titration with NaOH	99.0 - 101.0 %
GC (area %)	≥ 99.0 %

Specification: PRD.1.ZQ5.10000052620

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.