

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sorbitan monopalmitate

Product Number: 388920

MDL: MFCD00080946



TEST	Specification
Appearance (Color)	White to Yellow
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Carbon	59.1 - 72.2 %
Water (by Karl Fischer)	≤ 1.00 %
Solubility (Turbidity)	Clear
c = 50 mg/ml (5%); Ethanol	
Solubility (Color)	Colorless to Yellow

Specification: PRD.0.ZQ5.10000102012

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.