


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Methyl 1,3-dimethyl-1H-pyrazole-5-carboxylate

Product Number: **732664**
CAS Number: 10250-59-6
Formula: C₇H₁₀N₂O₂
Formula Weight: 154.17 g/mol



| TEST | Specification |
|------|---------------|
|------|---------------|

| | |
|--------------------|--|
| Appearance (Color) | Colorless to Yellow and Colorless to Brown |
| Appearance (Form) | Liquid |
| Infrared spectrum | Conforms to Structure |
| Purity (GC) | ≥ 88.0 % |

Specification: PRD.0.ZQ5.10000055870

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.