


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Acetylphenyl isocyanate - 97%

Product Number: 439932
CAS Number: 49647-20-3
MDL: MFCD00013879
Formula: C₉H₇NO₂
Formula Weight: 161.16 g/mol



TEST	Specification
Appearance (Color)	Yellow to Orange
Appearance (Form)	Solid or Liquid
Infrared spectrum	Conforms to Structure
Carbon	64.7 - 69.4 %
Nitrogen	8.3 - 9.0 %
Purity (GC)	≥ 96.5 %

Specification: PRD.0.ZQ5.10000049901

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.