


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Octanoic acid - $\geq 98\%$

Product Number: O3907
CAS Number: 124-07-2
MDL: MFCD00004429
Formula: C₈H₁₆O₂
Formula Weight: 144.21 g/mol



TEST	Specification
Appearance (Color)	Colorless or White
Appearance (Form)	Solid or Liquid
Infrared spectrum	Conforms to Structure
Purity (Titration by NaOH)	98.0 - 102.0 %
GC (area %)	$\geq 98.0\%$

Remarks:

Reviewing Chemist :

Specification: PRD.0.ZQ5.10000048021

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.