


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Potassium disulfate – ReagentPlus®, 99%

Product Number: 307750
CAS Number: 7790-62-7
MDL: MFCD00011385
Formula: K₂S₂O₇
Formula Weight: 254.32 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form) Powder and/or Chunks	Conforms to Requirements
Purity (Titration by NaOH)	≥ 98.5 %
ICP Major Analysis Confirms Potassium and Sulfur Components	Confirmed

Specification: PRD.0.ZQ5.10000032035

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.