

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Oncostatin M human - OSM, recombinant, expressed in HEK 293 cells, HumanKine®, suitable for cell culture

Product Number: H6541

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Purity by SDS-PAGE	≥ 95 %
Endotoxin Assay	< 1 EU/ug
Activity	0.1 - 1.5
ng/mL EC50	
The activity is determined	
by the dose-dependent stimulation of the proliferation	
of human TF-1 cells (human erythroleukemic indicator cell line)	

Specification: PRD.3.ZQ5.10000086159

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

