


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Poly(ethylene glycol) methyl ether amine - average  $M_n$  500

Product Number: **767565**  
CAS Number: 80506-64-5



Formula: CH<sub>3</sub>O(C<sub>2</sub>H<sub>4</sub>O)<sub>n</sub>C<sub>2</sub>H<sub>4</sub>NH<sub>2</sub>

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	Colorless to Dark Yellow and Colorless to Dark Orange
Appearance (Form)	Liquid
Proton NMR spectrum	Conforms to Structure
Molecular Number by HNMR	400 - 600
Poly Dispersity MW/MN	≤ 1.2

Specification: PRD.1.ZQ5.10000039400

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.