


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Amyloid β -Protein Fragment 35-25 - $\geq 95\%$ (HPLC)

Product Number: **A2201**
CAS Number: 147740-73-6
MDL: MFCD00214225
Formula: C₄₅H₈₁N₁₃O₁₄S
Formula Weight: 1,060.27 g/mol
Storage Temperature: -20 °C



TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
1 mg/mL in Acetonitrile:H ₂ O (7:3)	
% Peptide content	≥ 65
Amino acid analysis	
% Purity	≥ 95
by HPLC	

Specification: PRD.1.ZQ5.10000058676

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.