

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Disialoganglioside GD_{1a} from bovine brain - ≥95% (TLC), lyophilized powder

Product Number:

G2392

CAS Number:

12707-58-3

TEST	Specification
Appearance (Form)	Lyophilized Powder
Appearance (Colour)	White to Yellow/Tan
Solubility (Solvent)	CHCl ₃ /MeOH 2/1
Solubility (Conc)	19.60 - 20.40 mg/ml
Solubility (Heating)	Yes or No
Solubility (Turbidity)	Clear to Slightly Hazy
Solubility (Color)	Colorless to Light Yellow
Purity (TLC)	≥ 95.00 %
Recommended Retest Period	-----
3 Years	

Specification: PRE.0.ZQ5.10000001308

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.