

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Tris(2,2'-bipyridyl)dichlororuthenium(II) hexahydrate - powder

 Product Number:
 224758

 CAS Number:
 50525-27-4

 MDL:
 MFCD00149670

Formula: C30H24Cl2N6Ru · 6H2O

Formula Weight: 748.62 g/mol

TEST Specification

Red, Orange or Red-Orange

Appearance (Form)

Powder or Crystalline Powder

Infrared spectrum

Wave Length

Appearance (Color)

C = 0.01g/l in H2O

Extinction Coefficient

Carbon

Nitrogen

Solubility (Turbidity)

1% in H2O Solubility (Color)

Orange to Very Dark Orange, Red-Orange or Red

Specification: PRD.0.ZQ5.10000024693

Conforms to Requirements

Conforms to Requirements

Conforms to Structure

449 - 455 nm

> 11500

45.7 - 50.5 %

10.7 - 11.8 %

Clear to Slightly Hazy

Conforms