

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2-Fluoro-5-(trifluoromethyl)benzoic acid - 98%

Product Number:

455253

CAS Number:

115029-23-7

MDL:


MFCD00061292

Formula:

C₈H₄F₄O₂

Formula Weight:

208.11 g/mol



TEST

Specification

Appearance (Color)

White to Off-White

Appearance (Form)

Powder or Crystal or Chunk(s)

Infrared spectrum

Conforms to Structure

Purity (Titration by NaOH)

97.5 - 102.5 %

GC (area %)

≥ 97.5 %

Specification: PRD.0.ZQ5.10000023729

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.