

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification


Product Name:
Boc-3-iodo-D-Ala-OMe - 99%

Product Number: **426032**

MDL: MFCD00216580

Formula: C₉H₁₆INO₄

Formula Weight: 329.14 AMU



TEST

Specification

Appearance (Color)	White to Yellow
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Carbon	32.3 - 33.3 %
Nitrogen	3.9 - 4.5 %
Optical Rotation	3.0 - 5.0 deg
c = 2%; Methanol	

Specification Date : 07/06/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.