


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(1S,2R)-(+)-Norephedrine - 98%

Product Number: 317500
CAS Number: 37577-28-9
MDL: MFCD00064411
Formula: C₉H₁₃NO
Formula Weight: 151.21 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)	White to Light Yellow
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 97.5 %
Optical Rotation	38.0 - 42.0 deg
c = 7%; 1N HCl	

Specification: PRD.0.ZQ5.10000035927

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.