


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Xylometazoline hydrochloride

Product Number: **X6000**
CAS Number: 1218-35-5
MDL: MFCD00238707
Formula: C₁₆H₂₄N₂ · HCl
Formula Weight: 280.84 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
400 mg plus 4 mL of EtOH	
Purity (Titration by HCL04)	≥ 99 %
Purity (TLC)	≥ 99 %

Specification: PRD.0.ZQ5.10000049485

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.