

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

LB Broth with agar (Miller) - Highly-referenced nutrient-rich microbial growth powder medium with Agar, suitable for regular E.coli culture.

Product Number:

L3147

TEST	Specification
------	---------------

Appearance (Color)

Off White to Light Yellow to Light Brown

Appearance (Form)

Powder

Solubility (Color)

Almost Colorless to Yellow

Solubility (Turbidity)

Clear to Slightly Hazy

37 mg/mL, H₂O

Specification: PRD.2.ZQ5.1000008217

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.