

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Carbamyl-β-methylcholine chloride - ≥99% (TLC), crystalline

Product Number: C5259
CAS Number: 590-63-6
MDL: MFCD00055224
Formula: C7H17ClN2O2
Formula Weight: 196.68 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Powder
Carbon	42.1 - 43.4 %
Nitrogen	13.8 - 14.5 %
Purity (TLC)	≥ 99 %
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
100 mg/mL, H2O	

Specification: PRD.0.ZQ5.10000043780

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.