

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Stahl Aerobic Oxidation ABNO solution – 0.04M ABNO in Acetonitrile, Solution for Oxidation of Primary and Secondary Alcohols

Product Number:

796557



Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	Colorless to Dark Yellow and Colorless to Dark Orange
Appearance (Form)	Liquid
<sup>1</sup> H NMR Spectrum	Conforms to Structure
Component #1 2,2'-Bipyridyl	0.15 - 0.25 M
Component #2 ABNO	0.01 - 0.10 M
Component #3 1-Methylimidazole	0.35 - 0.45 M

Specification: PRD.2.ZQ5.10000076433

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

