

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)


Email USA: [techserv@sial.com](mailto:techserv@sial.com)

Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Poly(ethylene glycol) - average  $M_n$  950-1,050

Product Number: **P3515**  
 CAS Number: 25322-68-3  
 MDL: MFCD00081839  
 Formula:  $(C_2H_4O)_nH_2O$



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Conforms
Waxy Semisolid or Powder	
Solubility (Turbidity)	Clear
100mg/ml in H <sub>2</sub> O	
Solubility (Color)	Colorless
Brookfield Viscosity	≥ 20.0 cps
50% in H <sub>2</sub> O at 25 Degrees Celsius	
Recommended Retest Period	-----
4 Years	

Specification: PRD.2.ZQ5.10000018676

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.