

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com


Product Specification

Product Name:
Lithium diisopropylamide solution - 1.0 M in THF/hexanes

Product Number: 774766

Formula: C₆H₁₄LiN

Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Clarity)	Clear to Slightly Hazy
Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Liquid
Titration Sec-Butanol	0.95 - 1.10 M
Note	-----
Haziness may form upon storage but does not affect specification. If precipitation is present, then gently mix precipitate back into solution under an inert atmosphere.	
Recommended Retest Period	-----
9 months	

Specification: PRD.3.ZQ5.10000034008

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.