

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


Potassium carbonate sesquihydrate – Vetec™ reagent grade

Product Number:

V900082

CAS Number:

6381-79-9



Formula:

K₂CO₃ · 1.5H₂O

Formula Weight:

165.24 g/mol

TEST

Specification

Appearance (Colour)

White

Appearance (Form)

Conforms to Requirements

Crystals, Powder, Granules or Crystalline Powder and/or Chunks

X-Ray Diffraction

Conforms to Structure

Titration by HCL

98.5 - 101.0 %

Loss on Drying

14.0 - 16.5 %

Insoluble Matter

≤ 0.01 %

Specification: PRE.0.ZQ5.10000002536

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.