

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

3 α ,21-Dihydroxy-5 α -pregnan-20-one - $\geq 95\%$

Product Number:

P2016

CAS Number:


567-02-2

Formula:

C₂₁H₃₄O₃

Formula Weight:

334.49 g/mol



TEST	Specification
Appearance (Form)	Powder
Appearance (Colour)	White to White w/ Yellow Cast
Solubility (Solvent)	Chloroform
Solubility (Conc)	19.60 - 20.40 mg/ml
Solubility (Turbidity)	Clear
Solubility (Color)	Colorless
Water (by Karl Fischer)	
Purity (TLC)	$\geq 95.00\%$
NMR (Solvent)	CDCl ₃
Identity by NMR	Consistent
Recommended Retest Period	-----
3 Years	

Specification: PRE.2.ZQ5.10000021241

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

