

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Tetrachloroethylene - anhydrous, ≥99%

Product Number:

371696

CAS Number:

127-18-4

MDL:


MFCD00000834

Formula:

C₂Cl₄

Formula Weight:

165.83 g/mol



TEST

Specification

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

GC (area %)

≥ 99.0 %

Water (by Karl Fischer)

≤ 0.002 %

Note

For 100 mL unit size, water content may exceed above specification; not to exceed 0.005%.

Residue on Evaporation

≤ 0.0005 %

Specification: PRD.2.ZQ5.10000025786

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.