

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Gold - foil, thickness 0.025 mm, 99.99% trace metals basis

 Product Number:
 268461

 CAS Number:
 7440-57-5

 MDL:
 MFCD00003436

Formula: Au

Formula Weight: 196.97 g/mol

Au

TEST Specification

Appearance (Color)

Gold

Appearance (Form)

Foil

X-Ray Diffraction

Thickness 0.025 mm

Purity

99.99% Based On Trace Metals Analysis

Trace Metal Analysis

Specification: PRD.0.ZQ5.10000039908

Conforms

Conforms

Conforms to Structure

Conforms

Meets Requirements

 \leq 150.0 ppm

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.