

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Poly(styrene-co-allyl alcohol) - average  $M_w$  ~2,200, average  $M_n$  ~1,200, allyl alcohol 40 mol %

Product Number:

**468134**


CAS Number:

25119-62-4

MDL:

MFCD00134021

Formula:

C<sub>11</sub>H<sub>14</sub>O**TEST****Specification**

Appearance (Color)

Colorless to Faint Yellow

Appearance (Form)

Beads

Infrared spectrum

Conforms to Structure

Molecular Weight by GPC

2000 - 3200

Molecular Number by GPC

900 - 1500

Specification: PRD.0.ZQ5.10000034418

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.