


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Hydrocortisone - BioReagent, suitable for cell culture

Product Number: **H0888**
CAS Number: 50-23-7
MDL: MFCD00011654
Formula: C₂₁H₃₀O₅
Formula Weight: 362.46 g/mol



TEST

Specification

Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Very Faint Yellow
Solubility (Turbidity)	Clear
50 mg/ml, CHCl ₃ :EtOH (1:1)	
EmM	15.5 - 16.1
In MeOH	
Wavelength	242 - 243 nm
Purity (HPLC)	≥ 98 %
Cell Culture Test	Pass
Recommended Retest Period	-----
5 years	

Specification: PRD.3.ZQ5.10000006557

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.