


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
2-Bromothiophenol - 97%

Product Number: 242659  
CAS Number: 6320-02-1  
MDL: MFCD00004827  
Formula: C<sub>6</sub>H<sub>5</sub>BrS  
Formula Weight: 189.07 g/mol



### TEST

### Specification

Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Titration by Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub>	96.5 - 103.5 %
Purity (GC)	≥ 96.5 %

Specification: PRD.0.ZQ5.10000023291

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.