

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Medroxyprogesterone 17-acetate - $\geq 97\%$ (HPLC)

Product Number:

M1629

CAS Number:

71-58-9

MDL:


MFCD00010483

Formula:

C₂₄H₃₄O₄

Formula Weight:

386.52 g/mol

**TEST****Specification**

Appearance (Color)

White to Off-White

Appearance (Form)

Powder

EMM

15.5 - 16.3

Wavelength

240 nm

Purity (HPLC)

 $\geq 97\%$

Solubility (Color)

Colorless to Yellow

Solubility (Turbidity)

Clear

50 mg/ml, CHCl₃

Specification: PRD.0.ZQ5.10000042713

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.