

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Mercury – electronic grade, 99.9999% trace metals basis

Product Number: 294594
CAS Number: 7439-97-6
MDL: MFCD00011035
Formula: Hg
Formula Weight: 200.59 g/mol

Hg

TEST	Specification
Appearance (Color) Silver	Conforms
Appearance (Form) metallic liquid	Conforms
Trace Metal Analysis	≤ 2.000 ppm
Purity 99.9999% based on trace metals analysis	Conforms

Specification: PRD.1.ZQ5.10000068646

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.