

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


Potassium phenoxymethyltrifluoroborate - contains 20% KBr, 97%

Product Number:

710881

CAS Number:

1027642-30-3



Formula:

C7H7BF3KO

Formula Weight:

214.03 g/mol

TEST**Specification**

Appearance (Color)

White to Off-White

Appearance (Form)

Conforms to Requirements

Powder or Crystalline Powder

Proton NMR spectrum

Conforms to Structure

Fluorine NMR Spectra

Conforms to Structure

Titration by AgNO3

≤ 20.0 %

KBr

Carbon

37.9 - 40.7 %

Corrected for KBr Content

Specification: PRD.0.ZQ5.10000054408

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.