

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Titanium - wire, diam. 0.81 mm, 99.7% trace metals basis

Product Number: 267902
CAS Number: 7440-32-6
MDL: MFCD00011264
Formula: Ti
Formula Weight: 47.87 g/mol

Ti

TEST	Specification
------	---------------

Appearance (Color)	Silver or Grey
Appearance (Form)	Wire
ICP Major Analysis	Confirmed
Confirms Titanium Component	
Diameter	Conforms
0.81 mm	
Trace Metal Analysis	≤ 3500 ppm
Purity	Conforms
99.7% Based On Trace Metals Analysis	

Specification: PRD.0.ZQ5.10000035663

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.