

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


Hexane – Vetec™ reagent grade, anhydrous, ≥95%

Product Number:

V900805

CAS Number:

110-54-3



Formula:

C₆H₁₄

Formula Weight:

86.18 g/mol

TEST	Specification
------	---------------

Appearance (Colour)

Colorless

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 94.5 %

Residue on Evaporation

≤ 0.005 %

Water (by Karl Fischer)

≤ 0.005 %

Specification: PRE.0.ZQ5.10000004046

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.