

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Tungsten(VI) oxide – powder,  $\leq 25 \mu\text{m}$ ,  $\geq 99\%$  trace metals basis

**Product Number:** 232785  
**CAS Number:** 1314-35-8  
**MDL:** MFCD00011466  
**Formula:** O<sub>3</sub>W  
**Formula Weight:** 231.84 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Light Yellow to Dark Yellow and Light Green to Dark Green
Appearance (Form)	Powder
Loss on Ignition	$\leq 1.0 \%$
Particle Size	$\leq 25$ micron
ICP Major Analysis	Confirmed
Confirms Tungsten Component	
Purity	Conforms
$\geq 99\%$ Based On Trace Metals Analysis	
Trace Metal Analysis	$\leq 10000.0$ ppm

Specification: PRD.0.ZQ5.10000082054

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.