

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium chlorate - ReagentPlus®, ≥99%

Product Number: 244147
CAS Number: 7775-09-9
MDL: MFCD00003479
Formula: ClNaO3
Formula Weight: 106.44 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

White to Off-White or Pale Yellow

Appearance (Form)

Crystals or Crystals with Chunks

Infrared spectrum

Conforms to Structure

Titration by KMNO4

99.0 - 101.0 %

Specification: PRD.0.ZQ5.10000027193

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.