

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Pig IgG (whole molecule) antibody produced in rabbit - whole antiserum

Product Number:

P0916

Storage Temperature:

-20 °C

TEST**Specification**

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Very Light Brown to Very Dark Brown and Dark Brown-Red to Very Dark Brown-Red
Protein by Biuret	40 - 100 mg/ml
Antibody Concentration(Method)	Quantitative Precipitin Assay
Antibody Concentration	2.8 - 4.2 mg/ml
ELISA - Titer	1:40,000-1:80,000
Recommended Retest Period	-----
5 Years	

Specification: PRE.0.ZQ5.10000020796

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.