

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Fluorinert™ FC-40

Product Number: **F9755**  
CAS Number: 86508-42-1  
MDL: MFCD00131095

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Liquid
BOILING POINT	158.0 - 173.0 °C
Viscosity	1.6 - 2.8
units is cSt	
Density	1.845 - 1.895 g/ml
at 25 Degrees Celsius	
Vendor Information	Confirmed
Fluorinert is a trademark of 3M Company.	
3M is a trademark of 3M Company.	

Specification: PRD.2.ZQ5.10000004436

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.