

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Lithium iodate - 97%

Product Number: 443964
CAS Number: 13765-03-2
MDL: MFCD00016189
Formula: LiIO_3
Formula Weight: 181.84 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Powder or Solid or Crystals
Titration by $\text{Na}_2\text{S}_2\text{O}_3$ % IO_3	92.8 - 99.5 %
Titration	96.5 - 103.5 %
ICP Major Analysis Confirms Lithium Component	Confirmed

Specification: PRD.0.ZQ5.10000099134

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.