


3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name: Progesterone - ≥99%

**Product Number:** CAS Number: MDL: Formula: Formula Weight: P0130 57-83-0 MFCD00003658 C21H30O2 314.46 g/mol



## TEST

Appearance (Color) Appearance (Form) Solubility (Color) Solubility (Turbidity) 50 mg/ml, CHCl3 Infrared spectrum Purity (HPLC) Specific Rotation C = 2, dioxane, 25 deg C Recommended Retest Period 2 years

Specification: PRD.0.ZQ5.10000039425

## Specification

White to Off White Powder Colorless to Light Yellow Clear

Conforms to Structure > 99 % 175 - 183 °

-----

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.