


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium acrylate - 97%

Product Number: 408220
CAS Number: 7446-81-3
MDL: MFCD00067207
Formula: C₃H₃NaO₂
Formula Weight: 94.04 g/mol



TEST

Specification

Appearance (Color)	White
Appearance (Form)	Powder
Proton NMR spectrum	Conforms to Structure
Titration with HClO ₄	96.5 - 103.5 %
Solubility (Turbidity)	Clear
c = 5% In Water	
Solubility (Color)	Colorless to Very Faint Yellow

Specification: PRD.2.ZQ5.1000005159

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.